REHABI^{PVT®}

Penile Vacuum-Trainer

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Indication

Post-operative training of the cavernous body after radical prostatectomy and after surgical interventions in the minor pelvis. To improve erectile function and reduce potential incontinence after the surgery.

Training of the cavernous body in mild IPP (penile fibromatosis) to reduce the degree of the curvature. The device should be used over a longer period of time.

Product information

The **REHABI^{PVT}**® is a medical aid used in the rehabilitation and preservation of erectile function following prostatectomy and surgical interventions in the minor pelvis. The **REHABI^{PVT}**® generates negative pressure in the cylinder around the penis. This causes the cavernous body to fill up with blood, resulting in stiffening of the penis. These erections are of paramount importance for the rehabilitation and preservation of erectile function.¹⁻³ Oxygen-rich arterial blood flows into the cavernous body. This improves oxygen supply to the vessels and exercises the smooth muscles of the cavernous body, which in turn transforms the tissue of the cavernous body into inelastic connective tissue, thus reducing the length of the penis, the circumference of the penis and counteracting the loss of erectile function.

The **REHABI^{PVT}**® comes with the following parts:

- 1 cylinder, transparent
- 1 pump head, manual
- 1 adapter for the body
- Storage bag
- Lubricant Sylk® (50 ml)
- Instruction manual


How to use

Training of the cavernous body can commence immediately after successful wound healing and once the rehabilitation phase is started. Even if the device is only used at a much later point in time, the training of the cavernous body is still effective and helpful in treating impotence.


Training of the cavernous body can be done multiple times a day for several minutes.

There is no risk due to strong negative pressure as the training device creates only limited negative pressure. On the pump head there is a green "Air" button. By pressing this button you can aerate the cylinder and the vacuum will disappear.

1) Prepare all parts of the **REHABI^{PVT}**.

2) Place pump head onto the transparent cylinder
(Side of cylinder without sealing ring).


3) Place the adapter on the opposite, open side (with the sealing ring) of the cylinder.

4) Moisten adapter on the open side of the cylinder so that it seals well and can be fitted. Remove any pubic hairs.

5) To use the device, insert the penis into the open end of the cylinder. Negative pressure is created inside the cylinder by leverage. After about 20–25 pumping movements, a noticeable negative pressure should have built up in the cylinder. The cylinder must be pressed firmly onto the body for this purpose. The negative pressure increases the flow of blood into the cavernous body (corpus cavernosum) and the glans. This creates an artificial erection, which must be maintained for several minutes with the negative pressure. The penis remains in the cylinder during this time.

- 6) When you want to end the training of the cavernous body, allow air to enter the cylinder by pushing the "Air" button (venting valve). The negative pressure and the erection disappear.

Training of the cavernous body can be carried out multiple times a day as part of rehabilitation to improve the arterial blood supply.

Erectile tissue training with the **REHABI^{PVT}** can be performed over several months.

Erectile tissue training is only for rehabilitation purposes. The **REHABI^{PVT}** can be extended into a manual vacuum pump with constriction rings (for sexual intercourse). Request our „**REHABI^{PVT} extension to MES**“.

Cleaning

The **REHABI^{PVT}** can be easily cleaned with warm water and mild soap. To clean, remove the pump head from the cylinder. Clean the outside of the pump head with a damp cloth and soap. Do not use any potent cleaning agents.

You can store all parts in the enclosed bag.

Guarantee

There is a 2-year guarantee on the pump head. The guarantee does not apply in case of improper use, e.g. dropping the device.

In case of complaints, please indicate the serial number (found on the pump head).

Material

- Pump head: polycarbonate
- Cylinder: polycarbonate
- Adapter: TPE

Dimensions

REHABI^{PVT} complete Cylinder: 200 x 53 mm

Contraindication

Erectile tissue training should not be performed by persons with sickle cell anemia, leukemia, bone marrow tumours or other diseases associated with changes in blood coagulation and which lead to hemorrhage and prolonged erections.

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Tel. +49 (6105) 20 37 20
service@medintim.de
www.medintim.de

Studien / Studies / Arařtırmalar / Onderzoeken / Études :

- 1) Raina R, Agarwal A, Ausmundson S, Lakin M, Nandipati KC, Montague DK, Mansour D, Zippe CD. Early use of vacuum constriction device following radical prostatectomy facilitates early sexual activity and potentially earlier return of erectile function. International Journal of Impotence Research (2006) 18, 77-81.
- 2) Lehrfeld T, Lee DI. The role of vacuum erection devices in penile rehabilitation after radical prostatectomy. International Journal of Impotence Research (2009) 21, 158-164.
- 3) Sommer F. Mehr Potenz dank Schwellkörpertraining. Uro-News (2006) 1, 48-55.

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KESSEL medintim GmbH