




MEDintim
personal healthcare



Urethral *closure band*

(Nuremberg model)

The urethral closure band is a circular compression cuff which is easily applied to the penis. A soft compression pad prevents medicines injected into the urethra (liquids, gels, contrast means, etc.) from running out and similarly prevents the inadvertent loss of urine in mild to moderately severe urinary incontinence.

Advantages:

- the soft, padded, pleasant to wear cuff causes less pressure pain than comparable penis clamps
- the Velcro fastener simplifies adjustment
- the compression effect is easily adjusted by the patient
- light to wear and small in size the device provides good wearing comfort when running, sitting, walking
- Disposable

Material:

frontside: polyamid (100 %)
backside: cotton (100 %)
foam rubber: polyester (100 %)

Contraindications:

If you are allergic to cotton, polyamid or polyester, please contact your healthcare provider before use of this device.

Width: 2 cm

Contents: 3/10 closure bands

EAN 401327300078 9 (10 pieces)

EAN 410327300080 2 (3 pieces)

Art. no.: MED MHR

Made in Germany

 Medical Device



Kessel Marketing & Vertriebs GmbH
Kelsterbacher Straße 28
64546 Moerfelden-Walldorf, Germany

Tel. +49 (0) 6105 20 37 20
Fax +49 (0) 6105 45 59 01
Email: service@medintim.de



www.medintim.de