

Active LoopRing



Adjustable constriction ring
– reusable

Medintim


The Active LoopRing is a device, used to treat men with erectile dysfunction or impotence. The adjustable construction ring is intended to improve erections by increasing the fullness and firmness of the penis.

The ALR could be used as well with the Active³® Erection System or Manual Erection System (Vacuumtherapy), if the user wants to open the ring post climax and let sperm flow into the vagina.

Users also use the ALR for mild incontinence to prevent dripping. The ALR was originally developed to retain a medicine administered into the urethra and prevent it from leaking out.

Application:

The Active LoopRing consists of an adjustable loop made of silicone rubber that is placed around the base of the penis and tightened. The Active LoopRing may be tightened or loosened to fit. An O-ring and ball-locking unit allow for easy adjustment of the tension. Reusable device.



Material:

- Loop: TPE (Thermoplastic Elastomere)
- Ring/Ball: ABS plastic

Content:

- one Active LoopRing (Color: transparent)
- Instruction for use
- Box: 65 x 93 x 19 mm

Art. Nr.: IA ALR
UDI-DI 04013273001922

Available in pharmacies, specialist retailers
or directly from Medintim.

Hygiene product
Made in Germany

KESSEL medintim GmbH
Nordendstr. 82–84
64546 Mörfelden-Walldorf
Germany

Tel.: (06105) 20 37 20
Fax: (06105) 20 37 221
E-Mail: service@medintim.de
www.medintim.com