



Love
stressless

Material Bumpi® Rings:
100 % Thermoplastic Elastomer (TPE)

Material storage case:
Polypropylen (PP)

EAN 4013273000529

Made in Germany

Manufacturer:
KESSEL medintim GmbH
Nordendstr. 82 – 84
64546 Mörfelden-Walldorf
Germany
E-Mail: export@medintim.de

www.medintim.de

Bumpi® Rings

Soft for him.
Painless for her.

06.07.2023

Medintim



Women with endometriosis, prolapse of the uterus, after hysterectomy or with myomas sometimes experience bleeding or pain during sex. The medical name for this problem is collision dyspareunia.

Despite a variety of causes, there is a simple solution: the soft Bumpi® rings are slipped over the penis and in this way prevent the partner penetrating too deep. The result: pain-free sex, no more anxiety about injuries and finally relaxed, fulfilling sexuality.

Advantages:

- friendly to the skin, soft, flexible
- compatible with water-soluble lubricants, silicone oil, skin creams, vegetable fats
- biocompatible, tested to ISO 10993
- compatible with condoms made from latex and polyisoprene
- free from latex, phthalate, bisphenol A and silicone
- ecologically designed
- durable and reusable

Application:

One or several soft Bumpi® rings are stretched with both hands, and, before penetration, pulled up the penis shaft right up to the testicles.

The depth of penetration can be adjusted individually, depending on the number of Bumpi® rings used.

To pull off, the rings are again stretched and can then be removed.

Using the Bumpi® ring is simple, quick and does not interfere with lovemaking. There is no risk of the Bumpi® rings slipping inside the partner.

The Bumpi® rings lead to painless penetration and relaxed sexuality for both.

