


# Caya®

Contoured  
diaphragm



## Product benefits

- easy to insert and remove
- one size fits most women
- comfortable and practical for both partners

## Product specifications

- Elastic body: silicone rubber, biocompatible (according to ISO 10993)
- Core: contoured, one-piece polymer spring
- overall length: 75 mm
- overall width: 67 mm

## Product features

### Dynamic contours:

- folds in optimal insertion shape
- embeds gently in supportive tissue

### Gentle dynamics:

- soft and easy folding during insertion and removal
- slides gently out when removed

Each Caya® contoured diaphragm is individually checked for potential failures (100 % control). Produced and controlled in Germany.

### Manufacturer:

**KESSEL medintim GmbH**  
Nordendstr. 82 – 84  
64546 Mörfelden-Walldorf  
Germany

Tel.: +49 6105 20 37 20  
E-Mail: [service@medintim.de](mailto:service@medintim.de)

[www.medintim.de](http://www.medintim.de)

**CE** 0483  
Medical device  
Made in Germany

**Medintim**