

Caya® fits you, too.

The Caya® Diaphragm is a contraceptive for all women who want natural and safe contraception.

It is **especially** suitable for women who:

- want or need hormone-free contraception (e. g. in case of increased risk of thrombosis)
- cannot or do not want to use an IUD
- are breastfeeding
- want contraception only when they have sex
- have given birth / are in the interval between two pregnancies
- have a latex allergy
- want combination with natural methods of family planning (e. g. sympto-thermal method)

Contraindications:

- within the 6 weeks after childbirth
- in case of marked prolapse of pelvic floor / uterus
- if conventional Diaphragm sizes of 60 mm or 85 mm and more were used

Note:

Caya® does not protect against sexually transmitted diseases.

It is recommended to always use Caya® in combination with the Caya® gel.

Benefits of Caya® Diaphragm

- hormone-free
- easy to use
- no side effects
- reusable for two years
- latex-free
- no fitting required
- ecological
- no impairment of sexual feelings of both partners
- to be applied only when needed
- together with the diaphragm gel it provides **double safety:**

1. cervical cup as mechanical barrier
2. cellulose in diaphragm gel as additional mechanical barrier

Together with the Caya® Diaphragm you get:

- instruction for use
- storage case



Caya®

Contoured diaphragm



Caya®

Diaphragm gel



hormone-free
safe
easy to use



Caya®

Contoured diaphragm

Manufacturer:
KESSEL medintim GmbH
Kelsterbacher Str. 28
64546 Mörfelden-Walldorf
Germany
e-Mail: export@medintim.de
www.medintim.de

CE 0483
Medical device
Made in Germany

17.04.2023





Medintim

Contraception using the Caya® Diaphragm



The Caya® Diaphragm is a new contraceptive barrier device.

It is **anatomically** and **ergonomically** shaped and thus fits most women.



Many women wish for natural and safe contraception – without hormones!

The new Caya® Diaphragm offers reliable and self-determined contraception precisely when you need it. Caya® is inserted into the vagina before sex. It covers the cervix completely and is located between the posterior fornix and the niche on the pubic bone. During and after sex, Caya® prevents the sperm from entering the uterus. After intercourse, the Caya® Diaphragm must remain in the vagina for at least 6 hours.

For more information please see:
www.caya.eu



Double safety with Caya®

Caya® Diaphragm and Caya® Gel ensure safety:

1. The cervical cup is the mechanical barrier directly before the cervix. The solid rim reliably covers the cervix.
2. The Caya® gel contains cellulose, among other things, which forms an additional mechanical barrier directly at the cervix.

Safety, ease of use and wearing comfort of Caya® have been demonstrated in clinical studies.¹⁻⁴



References

1 Yang CC, Maravilla K, Kilbourne-Brook M, Austin G, et al. Magnetic resonance imaging of SILCS Diaphragm :anatomical consideration and corroboration with clinical fit. Contraception 2007; 76:238-244

2 Schwartz J, Ballagh S, Creinin M, et al. SILCS Diaphragm: postcoital testing of a new single-size contraceptive device. Contraception 2008; 78:237-244

3 Co ey P, Kilbourne-Brook M, Brache V, Cochón L, et al. Comparative acceptability of the SILCS and Ortho ALL-FLEX® Diaphragms among couples in the Dominican Republic. Contraception 2008; 78:418-423

4 Schwartz, Jill e.al 2015 Contraceptive Efficacy, Safety, Fit, and Acceptability of a single-size Diaphragm Developed with end-user input.