



Caya[®] pelvic model

Demonstration model of the female pelvis

An interactive teaching aid for illustrating the position of the **Caya[®]** contoured diaphragm in the female pelvis, for practising the insertion technique and palpation of the cervix.

This model is unique and was created to be anatomically correct. The organs of the small pelvis (bladder, uterus and rectum) are made of soft, highly elastic material. They can be removed from the pelvic floor after the transparent cover has been taken off.

The model separates into three parts.

Suitable for schools, family-planning centres, hospitals and gynaecology clinics.

The **Caya[®]** pelvic model can either stand freely, be laid flat or mounted on the wall.

Dimensions

Length 32.0 cm
Width 25.7 cm
Height 11.6 cm
Weight 1.5 kg

The **Caya[®]** pelvic model comprises the following components:

- Demonstration model of the female pelvis
- Instructions for use
- One **Caya[®]** contoured diaphragm (sample)

Cleaning

- Pelvic model can be easily cleaned with a damp cloth and mild soap

Caya[®] pelvic model
Art.-Nr. MO CPM
EAN 4013273001489

A solid storage case (**Caya[®]** pelvic model case) is also available for safekeeping and transportation.
Caya[®] pelvic model case (white)
Art.-Nr. IA KO W
EAN 4013273001496

© Original SOMSO[®] Model

Made in Germany