



MEDintim
personal healthcare



100 %
waterproof

ACTIVE *Sinnflut*[®] Man

Personal Massager for men – rechargeable

The **ACTIVE *Sinnflut*[®]** man is a new high technology concept for intimate relaxation of men in cases of light erection problems.

- waterproof
- extremely noiseless
- infinitely variable
- additional suction effect when covering the smaller hole
- soft surface due to medical grad silicone
- no batteries needed
- rechargeable due to induction
- strong motor and solid electronic

Instruction to use:

The **ACTIVE *Sinnflut*[®]** man device provokes and supports by its vibrating features a natural stimulated erection.

It helps as well in case of light erectile dysfunction.

Technical data:

material: ABS, Medical soft silicone
weight: 400 g
length: 205 mm
size: 82 x 65 mm
opening: 36 mm
deepness: 60 mm
colour: blue
loading station: ABS plastic with 220 V
power supply plug

Warranty:

two years warranty

Features:

On/off switch at the little button at the bottom of the device. With the little mini joystick you may increase or decrease or switch to impulse vibrations.

EAN: 4013273100014

Art.No.: VIB ASM

Made in Germany



KESSEL medintim GmbH
Nordendstr. 82-84
64546 Moerfelden-Walldorf, Germany

Tel. +49 (0) 6105 20 37 20
Fax (06105) 20 37 221
Email: service@medintim.de

www.medintim.de